

STERIS Applied Sterilization Technologies

The development in radiation sterilizing methods for single-use systems

ALFLOW Seminar 2023



Speaker

Pierre Reppert

Radiation Technology Centers Senior Manager EMEA-APAC

STERIS Applied Sterilization Technologies



STERIS Applied Sterilization Technologies

RADIATION MODALITIES



GAMMA

Exposes product to Cobalt 60 radiation



ELECTRON BEAM

Exposes product to high-energy electrons



X-RAY

Uses ionizing energy from electron beams

GAS MODALITIES



ETHYLENE OXIDE (EO)

Exposes product to gaseous sterilant



VAPORIZED HYDROGEN PEROXIDE (VHP)

Low temp gas process under deep vacuum

TESTING SERVICES



LABORATORY TESTING

Provides microbiological and analytical testing





PRODUCT & PACKAGE TESTING

Provides testing options for the validation of medical devices




Agenda

- Radiation Technology comparison
- Considerations
- Solutions: Risk based approach
- Publications


Technology Comparison: Gamma







Technology Comparison: Electron Beam



E-beam



Technology Comparison: X-ray



BPSA Publication




X-RAY STERILIZATION OF SINGLE-USE BIOPROCESS EQUIPMENT: PART I – INDUSTRY NEED, REQUIREMENTS AND RISK EVALUATION

2021

$$\Delta \frac{dH}{dt} = \Delta C_p \frac{dT}{dt} \quad E = h\nu = h \frac{c}{\lambda}$$

X-RAY STERILIZATION OF SINGLE-USE BIOPROCESS EQUIPMENT
PART I – INDUSTRY NEED, REQUIREMENTS AND RISK EVALUATION

BPSA
Bio-Process Systems Alliance
Advancing Single-Use Worldwide



Technology Comparison: Dose is Dose


Ionizing Radiation Sources - Dose given by secondary electrons

- All three technologies create ionization of materials
- Each ionization results in electrons distributed in product
- Electrons created ionize additional secondary ionizations further distributing through the product
- Regardless of ionizing source, electrons which are created inside the matters are doing the work

Gamma (photons)

X-Ray (photons)


E-beam (electrons)




→ **Dose in kGy**

Technology Comparison: Penetration


Photons (Gamma and X-ray)



Electron Beam



Technology Comparison: Penetration



Technology Comparison: Summary

Process comparison

Technology	Process	Process Load	Radiation Field Design
Gamma	<ul style="list-style-type: none">• Continuous• Batch• Incremental• Multi-pass	<ul style="list-style-type: none">• Tote• Pallet	<ul style="list-style-type: none">• Product overlapping• Source overlapping
X-ray	<ul style="list-style-type: none">• Continuous• Batch• Incremental• Multi-pass	<ul style="list-style-type: none">• Pallet	<ul style="list-style-type: none">• Product overlapping• Source overlapping
Electron Beam	<ul style="list-style-type: none">• Continuous• Batch	<ul style="list-style-type: none">• Box• Carrier	<ul style="list-style-type: none">• Source overlapping


Critical processing parameters:

- **Dose rate (kGy/h):** Risk of undesirable molecule recombination
- **Exposure time:** Risk of undesirable molecule recombination, ozone impact
- **Process capability:** Maximal dose received to achieve minimal required dose
- **Irradiation temperature:** Impact on some active material

Technology Comparison: Summary

	Gamma	E-beam	X-ray
Mode of action	Isotropic Photons	Electron	Photons with almost the same direction
DUR (~ 0.20 g/cc)	Dose range achievable : 25-40 kGy Ideal: 25-50 kGy	Dose range achievable : 25-50 kGy Ideal: 25-60 kGy	Dose range achievable : 25-35 kGy Ideal: 25-40 kGy
Dose rate	Variable during process A few kGy/h	Constant during process A few thousand kGy/h	Variable during process A few kGy/h to a few hundred of kGy/h
Temperature	Depends on design and Cobalt activity Typically, maximum temperature can go to 45 to 50°C	Depends on power Typically, maximum temperature can go to 50°C	Depends on power and design Typically, maximum temperature can go to 35°C to 40°C

Consideration: Transfer activities



Consideration: Minimal or Sterilization dose transfer

Transfer of verification dose or sterilisation dose (Section 8.4, 11137-1):

- Product not containing liquid water (dry product):
 - ✓ Transfer between facility operating the same radiation source is permitted without assessment
 - ✓ Transfer to different type of radiation source: assessment required
- Product containing liquid water:
 - ✓ Gamma to gamma is permitted
 - ✓ Two electron or X-ray operating under identical operating condition: permitted
 - ✓ Transfer to different type of radiation source: assessment required


Microbiological effectiveness of all 3 technologies has been demonstrated as equivalent therefore:

→ ISO 11137: Assessment consist of performing a Dose Audit exercise Other: Minimal dose used at current technology can be used (publication)

Publications : Minimal / Sterilization dose

Radiation Sterilization: Dose Is Dose

Joyce M. Hansen, Niki Fidopiastis, Trabue Bryans, Michelle Luebke and Terri Rymer
AAMI (2020)




Studies on the comparative effectiveness of X-rays, gamma rays and electron beams to inactivate microorganisms at different dose rates in industrial sterilization of medical devices.

Brian McEvoy, Ana Maksimovic, Daniel Howell, Pierre Reppert, Damien Ryan, Neil Rowan, Herve Michel (2023)
Radiation Physics and Chemistry

Microbicidal effectiveness of X-rays used for sterilization purposes

Tallentire, A. and Miller, A. (2015)
Radiation Physics and Chemistry



Consideration: Maximal acceptable dose transfer

ISO 11137-1: 8.1.1, The maximum acceptable dose for product shall be established. When treated with the maximum acceptable dose, product shall meet its specified functional requirements throughout its defined lifetime.

This statement is not specific to X-Ray, Gamma or EBeam. It applies to all three.

Clause 8.4.1 states that an assessment must be made to ensure that differences in conditions of two radiation modalities do not affect the validity of the maximum acceptable dose.

→ Assessment to be made → Risk Based Approach

Solutions: Risk based approach – Example 1

Critical processing parameters:

- **Dose rate:** Risk of undesirable molecule recombination
- **Exposure time:** Risk of undesirable molecule recombination, ozone impact
- **Process capability:** Maximal dose received to achieve minimal required dose
- **Irradiation temperature:** Impact on some active material

Processing Parameters	Gamma	X-ray	Electron beam
Dose rate	3	2	1
Exposure time	3	2	1
Maximal dose	2	1	3
Minimal dose	0	0	0
Irradiation temperature	3	1	2
Process capability	2	1	3
Total	13	7	10


0 Equivalent, 1 Best; 2 Medium; 3 worst

Solutions: Risk based approach – Example 2

Technology transfer - photons to photons Using a risk-based approach

	X-ray vs Gamma	Risk Evaluation
Sterilization / Minimal Dose	Equivalent	No
Maximal routine dose	Lower or Equivalent	No
Dose rate (kGy/h)	higher	Low
Temperature during process	Lower or Equivalent	Low
Exposure Time	Less	Low
Penetration (DUR)	Better or Equivalent	No
Product handling	Equivalent	No

Publications: AAMI Tir17: Compatibility of Materials subject to sterilization



Publications : Maximal Acceptable dose




X-ray: An effective Photon

Brian McEvoy, Hervé Michel, Daniel Howell and Philip Roxby

AAMI 2020

X-RAY STERILIZATION OF SINGLE USED BIOPROCESS EQUIPMENT Part I: Industry need, Requirements and Risk Evaluation

BPSA (2021)



Regulatory Approach for Transitioning from Gamma Ray to X-ray Radiation Sterilization

Alan Montgomery, Romain Bolle-Reddat, Shari Formica, Bradley Lundahl and Gerald McDonnell (AAMI 2021)

X-ray sterilization of biopharmaceutical manufacturing equipment

Roberto Menzel, Samuel Dorey, Tanja Maier, Ina Pahl, Armin Hauk

Biotechnology Progress 2021



Consideration: Activation assessment

Activation

ISO 11137-1 5.1.1 asks for an evaluation of a potential activation of materials with X-ray irradiation exceeding 5 MeV, or E-Beam with a >10 MeV E-beam treatment, even if the risk of activation of product is very small.

→ All products (Polymers, Implants, Animal feeds, API) tested to date at maximal acceptable dose have been declared as non-activated.

Tests were performed at SUVA (Swiss Government Accredited Laboratory) or at STERIS Libertyville RTC



Potential Induced Radioactivity in Material Process with X-ray Energy above 5 MeV

Hervé Michel, Thomas Kroc, Brian McEvoy, Deepak Patil, Pierre Reppert and Mark Smith – 06/2021

Consideration: Routine process transfer

Transference of Performance Qualification (mapping)

The transfer of performance qualification results is not possible as almost each installation and site are different

Statement valid for all 3 technologies

→ **PQ must be performed in accordance with ISO 11137 to assure the product specification can be met.**

Conclusions :

Critical processing parameters:

- **Dose rate (kGy/h):** Risk of undesirable molecule recombination
- **Exposure time:** Risk of undesirable molecule recombination, ozone impact
- **Process capability:** Maximal dose received to achieve minimal required dose
- **Irradiation temperature:** Impact on some active material

Considerations:

- Should the transfer from Gamma (photons) to X-ray (photons) be considered as a technology transfer?
- Can Gamma Qualification be considered as a worst-case scenario?

Solutions:

- Risk Based Approach
- Publications



Thank You!

Don't hesitate if you have questions in the future...
We are here to help!

Pierre Reppert

STERIS AST

RTC Senior Manager EMEA-APAC

Email: pierre_reppert@steris.com